

Choose your obturation solutions



ProTaper Gold®
paper points



Conform Fit™ gutta-percha
points for Protaper Gold®



Thermafil® for
ProTaper Gold®

INL / B EN PTGO TIP 000 / 03 / 2017 -
updated 11 / 2017



Mailefer Instruments Holding Sàrl
Chemin du Verger 3 - 1338 Ballaigues - Switzerland
Phone: +41 (0)21 843 92 92 - endo@dentsplysirona.com
www.dentsplysirona.com



**Dentsply
Sirona**
Endodontics




+
WE
KNOW
ENDO.


ProTaper Gold®

The same winning technique with **greater flexibility**

 **Dentsply
Sirona**
Endodontics

Treatment





Product line

ProGlider® Glide Path

PG



Variable Taper

016

2%

ProTaper Gold® Shaping Files

S1



018

2%

S2



020

4%

ProTaper Gold® Finishing Files

F1



020

7%

F2



025

8%

F3



030

9%

F4



040

6%

F5



050

5%

ProTaper Gold® SX (Accessory File)

SX



019

4%